

Discovery and characterization of C1qassociated mimotopes in SLE patients

Sarah Schott¹, Dagmar Hildebrand¹, Elke Hoffner¹, Fiordiligie Casilag¹, and Volker Stadler¹

BACKGROUND AND AIMS

A safe and effective SLE treatment should suppress the immune system as less as possible and target-oriented, to prevent side effects. Recently, antigen-specific treatment options came into focus. Peptides that mimic self-epitopes (mimotopes), block auto-antibodies and prevent mediated tissue damage are a promising tool to restore antigen-specific tolerance and to reduce clinical manifestations. Here we aimed to identify blocking mimotopes for anti-dsDNA autoantibodies and subsequently validate a possible harmful T-cell modulating effect of the mimotopes.



Figure 1. Graphical workflow. 1. For the identification of self-epitopes, sera are screened with PEPperCHIP® Peptide Microarrays. 2. Autoantibody-targeted epitopes are mimicked with synthetic peptides. Blocking values of peptides are evaluated in a competitive approach. 3. Functional blocking peptides are tested for mediated T-cell-activation.

METHODS

To identify self-epitopes, we designed a peptide microarray with non-nucleic acid targets of anti-dsDNA autoantibodies: 4309 linear peptides, 13 aa, overlap 10 aa, including 169 peptides of C1q. We analyzed the anti-IgG antibody response of SLE patients in comparison to healthy donors. Promising peptides were tested for specific blocking of antibodies and validated for T-cell activation via ELISpot assays (Fig.1).

¹PEPperPRINT GmbH, Heidelberg, Germany

RESULTS



Three of the C1q-derived peptides (IFDTVITNQEEPY and QGSEADSVFSGFL from the alpha subunit and GLPGLAGDHGEFG from the beta subunit) showed significantly more response in SLE patients and were selected for further validation (Fig.2).

FI <400 >400 <1000

QGSEADSVFSGFL

IFDTVITNQEEPY

GLPGLAGDHGEFG

2. Validation of blocking potential of C1q-associated peptides The mimotopes were synthesized and used in competitive assays. The sera of SLE patients were preincubated with the three peptides before the screen on the pep-

How can we support your research? Visit us at Get a free project consultation from our experts, or explore our website to learn more Booth 27 about our solutions for autoimmune research.

KEY FINDINGS

- Via PEPperCHIP® Peptide Microarrays we identified significant C1q-derived B-cell epitopes in SLE patients
- The epitope mimicking mimotopes blocked binding of serum-derived autoantibodies efficiently
- The mimotopes did neither induce cytotoxicity in PBMCs nor a harmful T-effector cell activation

1. Identification of B-cell self-epitopes in SLE patients

The analysis of the IgG responses revealed an increased expression of autoantibod-

tides with 0 values in all individuals were removed and 1795 peptides were further analyzed. (C) Heatmap of fluorescence intensities (FI) of 169 peptides of C1q. Black = FI below 500; yellow = F above 500; red = FI above 1000

tide microarray. All three mimotopes could effectively block binding of serum-derived autoantibodies to C1q-derived peptides.

3. Analysis of peptide-mediated T-cell activation The mimotopes were first tested for mediated cytotoxicity on PMCS. At this, no alarming change in cell viability was observed. Subsequently performed IFNg Fluorospots excluded a harmful activation of T effector cells (Fig.3).



Figure 3. (A) PBMCs were treated with mimotopres (50µg/ml) or staurossporine and stained with CTL-LDC[™] Live/Dead Cell Counting Kit. Viable cells: green; dead cells: red. (B) Results of an IFNg Fluorsopot assay. PBMCs were treated with solvent or pooled mimotopes or an CEF control Peptide Pool (10µg/ml).

SUMMARY AND CONCLUSION

We identified three C1q-derived self-epitopes in SLE patients. The respective soluble mimotopes blocked antibody-binding epitope-specific and showed no alarming activation of effector T-cells.

Peptide microarrays are a powerful tool to investigate the misdirected humoral immune response in autoimmune diseases and can support the development of peptide-based therapeutic strategies.

- PEPperPRINT GmbH Tullastr. 2 69126 Heidelberg Germany
- www.pepperprint.com ➡ info@pepperprint.com **L** +49-6221-7264489

in pepperprint-gmbh 🛛 🗙 @pepperchip 🕇 PEPperPRINTGmbH